Part no.	311-0099100-054
Product name	機器說明書/No Logo/英/.標準品/
Spec	L148*W574mm/雙面/模造紙/長邊彈簧6折 完成尺寸L148*W82mm/70P/黑
Designer	Dana
Color	C0 M0 Y0 K100
Info	





Snotty Boss logo here, like previous design

INTENDED USE

This product is intended for removal of the nasal secretions and mucus from children (0-12 years) at home environment. This nasal aspirator can be operated exclusively by adult. Children should not operate it.

IMPORTANT SAFETY PRECAUTIONS READ BEFORE USE

1. Use the device ONLY for the intended use described in this manual.

- 2. The device can remove nasal mucus from children 0-12 years (and
- physically compromised adults) in the comfort of home.
- 3. This product is **NOT** a toy and adult supervision and operation is required. 4. Keep the device and accessories away from young children. Some parts mav cause choking hazards.
- 5. Do **NOT** use accessories which are not included in this kit.
- 6. Do NOT use the device if faulty or damaged.
- 7. This device does **NOT** serve as a cure for any symptoms.
- 8. Do NOT use the device on other body parts.
- 9. Do **NOT** push the device too far into the nasal passage. OR use excessive force to insert the nozzle into the nose.
- 10.Do NOT disassemble or unscrew the device or attempt to repair it as factory warranty will be voided.
- 11.If your child's symptoms worsen consult a healthcare professional.
- 12.Do **NOT** use the device when patients have any of the following symptoms; nosebleed, acute head/facial/neck injury or coagulation/bleed-
- ing disorders. Please consult doctor if unsure. 13.To avoid eardrum pressure imbalance, during the use of the device, the
- children should not hold their breath, close their mouth or swallow. 14.Do **NOT** suction the child's nose continuously, but rather, off and on 10
- seconds at a time.
- 15. Please consult your healthcare professional about when and how often the nasal aspirator should be used.
- 16. If any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

KEEP THESE INSTRUCTIONS IN STORAGE BAG

CONTENTS

his kit contains: white handset 2 x AA alkaline batteries collection cup (2 pieces, top and base) 4 nozzles battery cover probe cleaning brush refillable saline spray bottle four-way probe for nose, ear, nail and navel cleaning storage bag

BEFORE YOU BEGIN

- Always use clean nozzles and collection cup to maintain hygiene and avoid cross-contamination
- Nozzles, collection cup can be removed from the device and washed after each use.
- If the patient feels any discomfort or any abnormality occurs during the process, stop aspiration immediately.
- Having spare nozzles is highly recommended in case of loss or damage.
- If the device is dropped or hit, check the device for damages or cracks, as it may affect function.

SPECIFICATIONS

- Model No.: TD-7601
- Power Source: Two 1.5V AA alkaline batteries
- Dimension & Weight: 41 (L) x 41 (W) x 225 (H) mm. 200 g
- Vacuum Pressure: 55 Kpa ~ 70 Kpa
- Noise Level: 45 dBA
- Ingress Protection: IP22
- Atmospheric Pressure Range: 700 hPA to 1060 hPA
- Operating Condition: 5°C to 40°C, 15% to 93% R.H.
- Storage Conditions: -25°C to 70°C, 10% to 95% R.H.

WARRANTY CLAIM

- You are entitled to a replacement or refund if the TD-7601 does not perform as described. This limited 12 month warranty applies to the TD-7601 only for motor malfunction unrelated to cleaning. Unscrewing any part of the aspirator will reconfigure the motor cog, and void warranty.
- The 60 day satisfaction guarantee applies if function does not meet expectations for nasal aspiration (full details on website).
- For faulty claims, please contact place of purchase, with receipt, for refund or exchange

QUICK HELP GUIDE

Scan QR code on unit base for easy troubleshooting and FAC





2. Insert two high quality 1.5V AA alkaline batteries to match the correct polarities (+ and -). Do NOT use rechargeshie batteries. 3. Firmly press the batteries to make sure they are completely inserted. 4. To place the cover back, just align the sides and snap it back.

- Always stand unit upright when full to avoid leakage into motor
- If snot flows into the blowhole/motor opening of the device body during use, the motor needs to be flushed, otherwise the motor may JAM or
- We recommend occasional internal motor flushing if snot leaks from

Only recommended if snot has leaked from Collection Cup or is visible

- Step 2: Power on the aspirator and dip the nozzle end into the basin. The aspirator will start to absorb water and then drain through the drain hole. Do not submerge whole aspirator body in water.
- Step 3: Flush for 90 seconds then remove nozzle from water and keep motor

Unflushed snot in the internal body may cause the motor cog movement to be affected by a dried piece of snot or cause internal rusting. Occasional internal flushing is essential as shown in diagram, however, DO NOT **UNSCREW** the unit to access internal parts as suction function will cease

TROUBLE SHOOTING

We want the TD-7001 to be the best device you have ever used to remove nasal mucus, and help your little one breathe, feed and sleep better.

- 1. Ensure batteries are new and inserted the right way.
- 2. Dry snot is impossible to suction, so squirt saline in the nostril first.
- 3. Ensure o-rings are in place and nozzles are fitted correctly.
- 4. Test the suction on palm of your hand. Remove top of Collection Cup. Press power button. Place palm of hand fully over circle opening on lower Collection Cup. You should feel good suction on your hand.
- Use a cool mist vaporiser in your child's room to hydrate and purify the air, all year round.

Consult your health professional if symptoms worsen or persist. Storage T. Scan QR code on unit base for easy troubleshooting and FAQ videos.

SYMBOL INFORMATION

- Storage conditions: room temperature in the enclosed carry bag.
- Avoid dropping or heavy impact.
- Avoid direct sunlight and high humidity.
- Remove batteries during extended storage.

Symbol	Referent	Symbol	Referent			
Ĩ	Consult instructions for use	Ŕ	Type BF applied part			
-	Manufacturer	X	Temperature limit			
SN	Serial number	Ì	Humidity limitation			
	Caution	IP22	Resistant to liquid ingress			
CE	CE mark	V _{RoHS}	RoHS Compliance			
UDI	Unique device identifier	MD	Medical device			
M	Manufacture date	*	Importer			
EC REP	Authorized representative in	n the Euro	pean Community			
X	This device does not belong to household waste and must be returned to a collection point for recycling electric and electronic devices according to local laws. If it contains batteries, the batteries should be removed and disposed in accordance with local regulations for separate collection of spent batteries.					

ECTREP MedNet EC-REP GmbH	Distributed by Snotty Noses Aust
Borkstraβe 10, 48163 Münster, Germany	Brisbane
TaiDoc Technology Corporation B1-7F, No. 127, Wugong 2nd Rd., Wugu Dist.,	snottynoses.com.a hello@snottynoses
24888 New Taipei City, Taiwan www.taidoc.com	CE A

stributed by	
notty Noses Australia Pty Ltd,	
isbane	
ottynoses.com.au •	
ello@snottynoses.com.au	
<u> </u>	
Ver. 5.0 2022/04	
23 MD / 311-0099100-054	

Manufacturer's declaration-electromagnetic immunity								
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.								
Immunity test	Immunity test IEC 60601 Compliance level test level							
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ±8 kV Air ±2 kV, ±4 kV, ±8 kV, ±15 kV	Contact: ±8 kV Air ±2 kV, ±4 kV, ±8 kV, ±15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.					
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input / output lines	Not applicable Not applicable	Mains power quality should be that of a typical home and professional healthcare environment.					
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	Not applicable Not applicable	Mains power quality should be that of a typical home and professional healthcare environment.					
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage dips: 0 % UT; 0,5 cycle 0 % UT; 1 cycle 70 % UT; 25/30 cycles Voltage interruptions: 0 % UT; 250/300 cycle	Voltage dips: Not applicable Not applicable Not applicable Voltage interruptions: Not applicable	Mains power quality should be that of a typical home and professional healthcare environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.					
Power frequency (50, 60 Hz) magnetic field IEC 61000-4-8 U	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz and 60 Hz	The device power frequency magnetic fields should be at levels characteristic of a typical location in a typical home and professional healthcare environment.					
NOTE UT is the a.c. m	ains voltage prior to appli	cation of the test level.						

The device is inte	ended for use in the election the user of the device sho	romagnetic environm	ent specified below.
The customer or		ould assure that it is u	used in such and environment.
Immunity test	IEC 60601	Compliance	Electromagnetic
	test level	level	environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms: 0,15 MHz – 80 MHz 6 Vrms: in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	Not applicable Not applicable	Portable and mobile RF communications equipment should be used no closer to a part of the device including cables, than the recommende separation distance calculate from the equation applicable the frequency of the transmit
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 800MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800MHz to 2,7 GHz Where P is the maximum out power rating of the transmitt watts (W) according to the transmitter manufacturer and the recommended separation distance in metres (m). Interference may occur in the vicinity of equipment marked with the following symbol: (b)

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Manufacturer's declaration-electromagnetic emission

The device is intended for use in the electromagnetic environment specified The customer or the user of the device should assure that it is used in such a

Emission test	Compliance	Electromagnetic er
RF emissions CISPR 11	Group 1	The device uses RF ene function. Therefore, its low and are not likely t in nearby electronic ec
RF emissions CISPR 11	Class B	The device is suitable establishments, includ
Harmonic emissions IEC 61000-3-2	Not applicable	to the public low-volta network that supplies
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	domestic purposes.

Recommended separation distance between portable and mobile RF communications equipment and tl

The device is intended for use in an electromagnetic environment in which ra disturbances are controlled. The customer or the user of the device can help electromagnetic interference by maintaining a minimum distance between communications equipment (transmitters) and the device as recommended the maximum output power of the communications equipment.

Rated maximum output power of	Separation distance according to freque (m)					
transmitter (W)	150 kHz to 80 MHz d =1,2√P	80 MHz to 800 MHz d =1,2√P				
0,01	N/A	0,12				
0,1	N/A	0,38				
1	N/A	1,2				
10	N/A	3,8				
100	N/A	12				

For transmitters rated at a maximum output power not listed above, the reco distance d in meters (m) can be estimated using the equation applicable to ti transmitter, where p is the maximum output power rating of the transmitter to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequ NOTE2 These guidelines may not apply in all situations. Electromagnetic pro absorption and reflection from structures, objects and people.

Normal Particular Particular<			lest specifications for E	NCLOSURE PORT IMMUN	IIY to RF wir	reless communicatio	ns equipment			Thank you for putting your trust in the 1D-7601. This booklet
understand servery method servery	The device is inte The customer or t	nded for use in the elect he user of the device sh	tromagnetic environment specified bel ould assure that it is used in such an en	ow. vironment.						provides important information to help you use, clean and maintain
energy of the server of the	Test frequency (MHz)	Band a)	Service ^{a)}	Modulation)	Maximum power (W)	Distance (m)	IMMUNITY TEST	Compliance LEVEL (V/m)	purchase docket/record your transaction details, for future
Image: Section of the sectin of the section of the section		(WITZ)					(,	(V/m)	(for home healthcare)	
d 459 49 <	385	380 - 390	TETRA 400	Pulse modulation ^{b)} 18 Hz		1,8	0,3	27	27	Data of susshares
100 74-787 LTE Band 13, 17 Note modulation* 217 Hz 0.2 0.3 9 9 100 800-960 6504 800/000, TETRA 800, DED 820, COM 8100, UTE Band 5 Note modulation* 18 Hz 2 0.3 28 28 1720 1700-1790 CSM 800/000, TETRA 800, UTE Band 5 Note modulation* 18 Hz 2 0.3 28 28 1720 1700-1790 CSM 1900, COM 8100, UTE Band 13, 4, 25 (UMTS) Pale modulation* 217 Hz 2 0.3 28 28 1845 115 Band 13, 4, 25 (UMTS) Pale modulation* 217 Hz 2 0.3 28 28 1940 2400-2570 Bluetooth, WLAN, 802.11 a/m Pale modulation* 217 Hz 0.2 0.3 9 9 1950 5000 WLAN 802.11 a/m Pale modulation* 217 Hz 0.2 0.3 9 9 1000-45. 5000 WLAN 802.11 a/m Pale modulation* 217 Hz 0.2 0.3 9 9 1000-45. 167 terme service, only the uplink frequencies are included. 10.1 Hz	450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz deviation 1 kHz sine		2	0,3	28	28	Place of purchase:
745 170 170 170 170 170 170 190 00 600 000 900 000 900 000 900 000 900 000 900 100 900 900 100 900 900 100 900 900 100 900 900 100 900 900 100 90	710	704 – 787	LTE Band 13, 17	Pulse modulation b		0,2	0,3	9	9	Receipt number/Order number:
780 Image: Constraint of the second seco	745			217 HZ						· · · · · · · · · · · · · · · · · · ·
810 800-900 CSM 800/900, TER 8400, UE Band S. Pule modulation* 2 0.3 28 28 28 1720 1700-1990 CSM 1000; CSM 1800; UE Band 1, 3, 4, 2; UMTS Pule modulation* 2 0.3 28 28 1845 1720 1700-1990 CSM 100; CSM 1900; UE Band 1, 3, 4, 2; UMTS Pule modulation* 2 0.3 28 28 1970 0 0 CSM 800; DECT; UE Band 1, 3, 4, 2; UMTS Pule modulation* 2 0.3 28 28 1970 0 100-500 WLAN 802; I1 br/m, BFD 250; UT Band 7 Pule modulation* 2 0.3 28 28 5500 100-500 WLAN 802; I1 br/m, BFD 250; UT Band 7 Pule modulation* 0.2 0.3 9 9 5785 100-500 WLAN 802; I1 br/m, 27 fz Pule modulation, T ME 07 THE STIEM may be reduced to 1 m. The 1 m test distance is permitted by EC Enotive Boss Enotive Boss 1000-51. 100-500; A an altenative to FM modulation, 50% duty cycle square wave signal. 1 Number of the start and and safety requirements of: IEC/EN 60601-1-1; IEC/EN 60601-1-1; I	780									
870 ITE Band 5 ITE Band 7	810	800 - 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850	Pulse modulation ^{b)}		2	0,3	28	28	
120 1700 - 1900 GM 1800 (CDM 100); Pule modulation * 2 0.3 28 28 1120 1700 - 1900 GM 1800 (CDM 400); Pule modulation * 2 0.3 28 28 1970 1970 Bluetooth, WLAN, 802.11 big/n, Pule modulation * 2 0.3 28 28 2450 2400 - 2570 Bluetooth, WLAN, 802.11 big/n, Pule modulation * 2 0.3 28 28 5500 17785 177 Hz 177 Hz 0.2 0.3 9 9 NOTE if necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC Broitiv Boss 1 10 or some services, only the uplik frequencies are included. Bit modulation, 30 % pulse modulation, 1 twould be worst case. This device has been tested to meet the electrical and safety requirements of: IEC//EN 60601-1.1 IEC//EN 60601-1.21. IEC 606001-1.21. Frequirement for safety IEC60001-1.2 Requirements for EMC 1 <td>870</td> <td>_</td> <td>LTE Band 5</td> <td>10112</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>	870	_	LTE Band 5	10112						
1720 1700-1990 GSM 1800; CDCM 1900; GSM 1800; CDCM Pulse modulation* 2 0.3 28 28 28 1970 2400-2570 Bluetooth, WLAN, 802.11 b/g/n RFID 2450, UTE Band 7 Pulse modulation* 2 0.3 28 28 28 5240 5100-5800 WLAN 802.11 a/n Pulse modulation* 2 0.3 28 28 28 5500 5500 WLAN 802.11 a/n Pulse modulation* 2 0.3 9 9 9 NOTE If necessary to achieve the MMUNITY TEST LEVEL, the distance between the transmitting antenna af the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 6:000-4.3. 3:05 % 0.007 cycle square wave signal. c) As an alternative to FM modulation, at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. This device has been tested to meet the electrical and safety requirements of: IEC/EM 60601-1, IEC/EM 60601-1-21. IEC 60601-1-11. Reference to Standards: • IEC60001-1-2. IEC 60601-1-11. Reference to Standards: • IEC60001-1-2. Requirements for EMC	930									Diligent external cleaning and internal motor flushing is
1845 ITE Band 1, 3, 4, 2; UMTS Image: Comparison of the com	1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT;	Pulse modulation ^{b)} 217 Hz		2	0,3	28	28	essential to keep your Tu-rout under warranty and in working
1970 1970 2450 2400-2570 Bluetooth, WLAN, 802.11 b/g/n, Pulse modulation ⁴⁴ 217 Hz 200 5240 5240 5200 5765 NOTE If necesary to achieve the MMUNITY TEST LEVEL, the distance between the transmitting antenna ard the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50 % duty cycle square wave signal. c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1-11. Reference to Standards: ·IEC606061-1-11. Standards: ·IEC606061-1-2. IEC0011-2. Requirements for safety ·IEC606061-1-2. IEC0011-2. Requirements for safety ·IEC00061-1-11. Standards: ·IEC00601-1-11. ·IEC00601-1-11. Reference to Standards: ·IEC00601-1-11. ·IEC00601-1-11. Reference to Standards: ·IEC006001-1-2. ·IEC00601-1-11.	1845	_	LTE Band 1, 3, 4, 25; UMTS							order
2450 2400 - 2570 Bluetooth, WLAN, 802.11 bg/n, Pulse modulation ^M 2 0.3 28 28 5240 5100 - 5800 WLAN 802.11 a/n Pulse modulation ^M 0.2 0.3 9 9 5785 NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 1000 -4.3. a) for some services, only the uplink frequencies are included. b) The carrier shall be modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1-2, IEC 60601-1-	1970									Bnotty Boss
5240 5100-5800 WLAN 802.11 a/n Pulse modulation * 0.2 0.3 9 9 5705 5785 NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna ard the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna ard the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC b) The carrier shall be modulated using a 50 % duty cycle square wave signal. c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-1-11. Reference to Standards: • IECC60001-1, IEC/EN 60601-1-2, IEC 60601-1-11. Reference to Standards: • IECC60001-1, IEC 60601-1-12, Requirements for safety · IECC60001-1, IEC 60601-1-12. Requirements for safety	2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz		2	0,3	28	28	
5500 5785 NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna ar d the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3. a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50% duty cycle square wave signal. c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1, IEC/EN 60601-1, IEC 60001-1, IEC 6000-1, IEC	5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz		0,2	0,3	9	9	
5785 NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna ar d the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 6100-4-3. a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50 % duty cycle square wave signal. c) A an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-1-11. Reference to Standards: • IEC608001-1.2 Requirement for safety • IEC60801-1.2 Requirement for EMC	5500	_		217112						
NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3. a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50 % duty cycle square wave signal. c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-1-11. Reference to Standards: • IEC600601-1, IEC 60601-1-12 Requirement for safety • IEC600601-1-2 Requirements for EMC	5785									
a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50 % duty cycle square wave signal. c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1-2, IEC 60601-1-11. Reference to Standards: · IEC60601-1. IEC 60601-1-2 Requirement for safety · IEC60601-1.2 Requirements for EMC	NOTE If necessary 61000-4-3.	to achieve the IMMUN	TY TEST LEVEL, the distance between the	ne transmitting antenna ar	d the ME EQU	JIPMENT or ME SYSTE	M may be reduced	to 1 m. The 1 m test dis	stance is permitted by IEC	
b) The carrier shall be modulated using a 50% duty cycle square wave signal. c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1-11. Reference to Standards: • IEC60601-1-1 General requirement for safety • IEC60601-1-2 Requirements for EMC	a) For some servic	es only the unlink freq	iencies are included							
This device has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-1-11. Reference to Standards: • IEC60601-1, IEC 60601-1-1 Requirement for safety • IEC60601-1, 2 Requirements for EMC	b) The carrier shal	be modulated using a	50 % duty cycle square wave signal.	od hosouso while it doos	at represent	actual madulation it				
hg he has been tested to meet the electrical and safety requirements of: IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-1-11. Reference to Standards: • IEC60601-1, IEC 60601-1-11 General requirement for safety • IEC60601-1-2 Requirements for EMC		e to FM modulation, 50	% puise modulation at 18 Hz may be u	sed because while it does in	iot represent	actual modulation, it	would be worst ca	se.		
By Reference to Standards: • IEC60601-1, IEC 60601-1-11 General requirement for safety • IEC60601-1-2 Requirements for EMC					This de IEC/EN	evice has been te I 60601-1, IEC/E	ested to meet 1 N 60601-1-2,	he electrical and s IEC 60601-1-11.	safety requirements of:	
					Referer • IEC60 • IEC60	nce to Standards)601-1, IEC 606)601-1-2 Require	s: 01-1-11 Gene ements for EM	ral requirement fo C	r safety	